

Fixture thread options for better initial stability

The option of different thread depths (regular or deep) and special KnifeThread® design enable easy implant placement with good primary stability in ALL bone densities


- REGULAR thread recommended for hard bone (D1 & D2)
- DEEP thread recommended for soft bone or poor bone density (D3 & D4)

	Fixture Diameter				
	Ø3.3	Ø3.7	Ø4.1	Ø4.4	Ø4.8
Regular Thread					
Thread Depth	0.4	0.4	0.45	0.45	0.4
Deep Thread					
Thread Depth	0.6	0.6	0.65	0.6	0.65



Regular & Deep Threads

Fixture Diameter	Apical Diameter (Regular & Deep Thread)	Core Diameter	Thread Depth (Regular Thread / Deep Thread)		Length(mm)	Connection Diameter
Ø3.3	Ø1.1	Ø2.8	0.4	0.6	7 / 7.7 / 9.2 / 10.7 12.2 / 14.2 / 17.2	Ø2.8
Ø3.7	Ø1.4	Ø3.2	0.4	0.6		Ø2.8
Ø4.1	Ø1.9	Ø3.5	0.45	0.65		Ø3.3
Ø4.4	Ø2.1	Ø3.8	0.45	0.6		Ø3.3
Ø4.8	Ø2.4	Ø4.2	0.4	0.65		Ø3.3



Fixture Sizes & Packaging

I. Fixture Sizes

NC Ø3.3 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø3.3	NC	7	ARO3307C
		8.5	ARO3308C
		10	ARO3310C
		11.5	ARO3311C
		15	ARO3315C
		18	ARO3318C



NC Ø3.7 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø3.7	NC	7	ARO3707C
		8.5	ARO3708C
		10	ARO3710C
		11.5	ARO3711C
		13	ARO3713C
		15	ARO3715C
18	ARO3718C		



RC Ø4.1 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.1	RC	7	ARO4107C
		8.5	ARO4108C
		10	ARO4110C
		11.5	ARO4111C
		13	ARO4113C
		15	ARO4115C
18	ARO4118C		



RC Ø4.4 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.4	RC	7	ARO4407C
		8.5	ARO4408C
		10	ARO4410C
		11.5	ARO4411C
		13	ARO4413C
		15	ARO4415C
18	ARO4418C		



RC Ø4.8 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.8	RC	7	ARO4807C
		8.5	ARO4808C
		10	ARO4810C
		11.5	ARO4811C
		13	ARO4813C
		15	ARO4815C
18	ARO4818C		



RC Ø5.3 Regular Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø5.3	RC	7	ARO5307C
		8.5	ARO5308C
		10	ARO5310C
		11.5	ARO5311C
		13	ARO5313C
		15	ARO5315C
		18	ARO5318C



NC Ø3.3 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø3.3	NC	7	ARO3307DC
		8.5	ARO3308DC
		10	ARO3310DC
		11.5	ARO3311DC
		13	ARO3313DC
		15	ARO3315DC
		18	ARO3318DC



NC Ø3.7 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø3.7	NC	7	ARO3707DC
		8.5	ARO3708DC
		10	ARO3710DC
		11.5	ARO3711DC
		13	ARO3713DC
		15	ARO3715DC
		18	ARO3718DC



RC Ø4.1 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.1	RC	7	ARO4107DC
		8.5	ARO4108DC
		10	ARO4110DC
		11.5	ARO4111DC
		13	ARO4113DC
		15	ARO4115DC
		18	ARO4118DC



RC Ø4.4 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.4	RC	7	ARO4407DC
		8.5	ARO4408DC
		10	ARO4410DC
		11.5	ARO4411DC
		13	ARO4413DC
		15	ARO4415DC
		18	ARO4418DC



RC Ø4.8 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø4.8	RC	7	ARO4807DC
		8.5	ARO4808DC
		10	ARO4810DC
		11.5	ARO4811DC
		13	ARO4813DC
		15	ARO4815DC
		18	ARO4818DC



RC Ø5.3 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø5.3	RC	7	ARO5307DC
		8.5	ARO5308DC
		10	ARO5310DC
		11.5	ARO5311DC
		13	ARO5313DC
		15	ARO5315DC
		18	ARO5318DC



RC Ø5.8 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø5.8	RC	7	ARO5807DC
		8.5	ARO5808DC
		10	ARO5810DC
		11.5	ARO5811DC
		13	ARO5813DC
		15	ARO5815DC
		18	ARO5818DC



RC Ø6.3 Deep Thread

- Cover Screw included

Fixture Diameter	Connection	Length (mm)	Ref.C
Ø6.3	RC	7	ARO6307DC
		8.5	ARO6308DC
		10	ARO6310DC
		11.5	ARO6311DC
		13	ARO6313DC
		15	ARO6315DC
		18	ARO6318DC

